

Evaluation of correlation between Diopsys® NOVA™ fixed-luminance flicker ERG and Diagnosys® Espion 2™ flicker ERG parameters

Jonathan Regenold^{1,2}, Hien Luong Doan¹, Hashem Ghoraba^{1,3}, Hassan Khojasteh^{1,4}, Jaclyn Joyce Hwang¹, Negin Yavari¹, Amir Akhavanrezayat¹, Ngoc Trong Tuong Than¹, Anthony Huy Dinh Le¹, Muhammad Sohail Halim¹, Quan Dong Nguyen¹

1. Spencer Center for Vision Research, Byers Eye Institute, Stanford University School of Medicine, 2452 Watson Court, Suite 200, Palo Alto, CA, 94303, USA.
2. University of Cincinnati College of Medicine, Cincinnati, OH, USA.
3. Cole Eye Institute, Cleveland Clinic, Cleveland, OH, USA.
4. Bascom Palmer Eye Institute, Miami, FL, USA.

Purpose: Diopsys® NOVA™ is a novel full-field electroretinography (ffERG) device that can make rapid measurements of retinal electrophysiologic function. Diagnosys® Espion 2™ is a clinical gold-standard ERG device. This study aimed to investigate whether light-adapted Diopsys® NOVA™ fixed-luminance flicker ffERG magnitude and implicit time (converted from phase) measurements correlate with light-adapted Diagnosys® Espion 2™ flicker ffERG amplitude and implicit time measurements, respectively.

Methods: Twelve patients (22 eyes) with various retinal and uveitic diseases underwent light-adapted Diagnosys® Espion 2™ and Diopsys® NOVA™ fixed-luminance flicker testing. Diopsys® magnitude and implicit time (converted from phase) measurements were compared to Diagnosys® amplitude and implicit time measurements, and a Pearson correlation was used to evaluate any existing correlation. Groups were also compared using generalized estimating equations. Bland-Altman plots were utilized to determine agreement between the comparison groups.

Results: Age of patients ranged from 14 to 87 years. 58% (n = 7/12) of patients were female. A significant, positive correlation ($r = 0.880$, $P < 0.001$) was observed between magnitude (Diopsys®) and amplitude (Diagnosys®) measurements. Amplitude increases by 6.69 μV for each 1 μV increase in Magnitude (p-value < 0.001). A statistically significant, strong positive correlation was observed between Diopsys® implicit time measurements (converted from phase) and Diagnosys® implicit time measurements ($r = 0.814$, p-value < 0.001). For each 1 ms increase in Diopsys® implicit time, Diagnosys® implicit time increases by 1.13 ms (p-value < 0.001).

Conclusions: There is a statistically significant positive correlation between lightadapted Diopsys® NOVA™ fixed-luminance flicker amplitude and Diagnosys® flicker magnitude values. Additionally, there is a statistically significant positive correlation between Diopsys® NOVA™ fixed-luminance flicker implicit time (converted from phase) and Diagnosys® flicker implicit time values. These results imply that the Diopsys® NOVA™ module, which utilizes the nonstandard shortened International Society for Clinical Electrophysiology of Vision (ISCEV) ERG protocol, can produce reliable lightadapted flicker ffERG measurements.

Keywords: Diagnosys®; Diopsys®; Fixed luminance; Flicker electroretinography; Fullfield electroretinography.

Evaluation of correlation between Diopsys® NOVA™ fixed-luminance flicker ERG and Diagnosys® Espion 2™ flicker ERG parameters

Jonathan Regenold^{1,2}, Hien Luong Doan¹, Hashem Ghoraba^{1,3}, Hassan Khojasteh^{1,4}, Jaclyn Joyce Hwang¹, Negin Yavari¹, Amir Akhavanrezayat¹, Ngoc Trong Tuong Than¹, Anthony Huy Dinh Le¹, Muhammad Sohail Halim¹, Quan Dong Nguyen¹

1. Spencer Center for Vision Research, Byers Eye Institute, Stanford University School of Medicine, 2452 Watson Court, Suite 200, Palo Alto, CA, 94303, USA. 2. University of Cincinnati College of Medicine, Cincinnati, OH, USA. 3. Cole Eye Institute, Cleveland Clinic, Cleveland, OH, USA. 4. Bascom Palmer Eye Institute, Miami, FL, USA.

Introduction

Light-adapted flicker electroretinography (ERG) is an electrophysiologic test of the retina that evaluates cone and bipolar cell function.¹ Diopsys® NOVA is a novel ERG modality that can perform a variety of ERG tests including 30 Hz flicker ERG. Diagnosys® Espion 2 is a gold-standard clinical ERG device. Due to its user-friendly design and affordability, the Diopsys® NOVA device has the potential to improve the accessibility of ERG testing.

Purpose

This study aimed to investigate whether light-adapted Diopsys® NOVA™ fixed-luminance flicker ffERG magnitude and implicit time (converted from phase) measurements correlate with light-adapted Diagnosys® Espion 2™ flicker ffERG amplitude and implicit time measurements, respectively.


Methods

Twelve patients (22 eyes) with various retinal and uveitic diseases underwent light-adapted Diagnosys® Espion 2™ and Diopsys® NOVA™ fixed-luminance flicker testing. Diopsys® magnitude and implicit time (converted from phase) measurements were compared to Diagnosys® amplitude and implicit time measurements, and a Pearson correlation was used to evaluate any existing correlation. Groups were also compared using generalized estimating equations. Bland-Altman plots were utilized to determine agreement between the comparison groups.

Results

Age of patients ranged from 14 to 87 years. 58% (n = 7/12) of patients were female. A significant, positive correlation ($r = 0.880$, $P < 0.001$) was observed between magnitude (Diopsys®) and amplitude (Diagnosys®) measurements. Amplitude increases by $6.69 \mu\text{V}$ for each $1 \mu\text{V}$ increase in Magnitude (p-value < 0.001). A statistically significant, strong positive correlation was observed between Diopsys® implicit time measurements (converted from phase) and Diagnosys® implicit time measurements ($r = 0.814$, p-value < 0.001). For each 1 ms increase in Diopsys® implicit time, Diagnosys® implicit time increases by 1.13 ms (p-value < 0.001).

Fig. 1



Diopsys® ERG Lid Sensor Diagnosys® DTL Plus Electrode

Fig. 1. Comparison between electrodes utilized by the devices, showcasing the patient-friendly design of the Diopsys® device

Fig. 2




Fig. 2. Bland–Altman plot analyzing agreement between Diopsys® magnitude and Diagnosys® amplitude measurements. Similar results were shown for time measurements.

Conclusions

There is a statistically significant positive correlation between light-adapted Diopsys® NOVA™ fixed-luminance flicker amplitude and Diagnosys® flicker magnitude values. Additionally, there is a statistically significant positive correlation between Diopsys® NOVA™ fixed-luminance flicker implicit time (converted from phase) and Diagnosys® flicker implicit time values. These results imply that the Diopsys® NOVA™ module, which utilizes the nonstandard shortened International Society for Clinical Electrophysiology of Vision (ISCEV) ERG protocol, can produce reliable light-adapted flicker ffERG measurements.

References

1. Creel DJ (2019) Electroretinograms. *Handb Clin Neurol* 160:481–493.
2. https://diopsys.com/press_releases/advanced-erg-lid-sensor-revolutionizes-eye-care-testing/
3. Brouwer et al. Effects of DTL electrode position on the amplitude and implicit time of the electroretinogram. *Documenta Ophthalmologica*. 2019.